

For Research Use Only

TRIM59 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60824-3-PBS



Basic Information

Catalog Number: 60824-3-PBS	GenBank Accession Number: NM_173084	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 286827	CloneNo.: 1F3C3
Source: Mouse	UNIPROT ID: Q8IWR1	
Isotype: IgG1	Full Name: tripartite motif-containing 59	
Immunogen Catalog Number: AG29572	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51219-2, TRIM59 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60824-3-PBS. Detection antibody: 60824-4-PBS. Standard:Ag29572. Range: 0.781-100 ng/mL.