

For Research Use Only

TMEM176B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60825-2-PBS



Basic Information

Catalog Number: 60825-2-PBS	GenBank Accession Number: BC004358	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 28959	CloneNo.: 2A9G7
Source: Mouse	UNIPROT ID: Q3YBM2	
Isotype: IgG1	Full Name: transmembrane protein 176B	
Immunogen Catalog Number: AG13918	Calculated MW: 270 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51220-1, TMEM176B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60825-1-PBS. Detection antibody: 60825-2-PBS. Standard:Ag13918. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51220-2, TMEM176B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60825-3-PBS. Detection antibody: 60825-2-PBS. Standard:Ag13918. Range: 0.098-100 ng/mL