

For Research Use Only

FHIT Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60833-3-PBS



Basic Information

Catalog Number: 60833-3-PBS	GenBank Accession Number: BC032336	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2272	CloneNo.: 1F7G3
Source: Mouse	UNIPROT ID: P49789	
Isotype: IgG1	Full Name: fragile histidine triad gene	
Immunogen Catalog Number: AG14471	Calculated MW: 147 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51232-2, FHIT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60833-2-PBS. Detection antibody: 60833-3-PBS. Standard:Ag14471. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51232-3, FHIT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60833-4-PBS. Detection antibody: 60833-3-PBS. Standard:Ag14471. Range: 0.391-100 ng/mL.