

For Research Use Only

WBP2NL Monoclonal antibody, PBS Only (Detector)

Catalog Number:60839-2-PBS



Basic Information

| | | |
|---|--|--|
| Catalog Number: 60839-2-PBS | GenBank Accession Number: BC022546 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 164684 | CloneNo.: 1F6A3 |
| Source: Mouse | UNIPROT ID: Q6ICG8 | |
| Isotype: IgG1 | Full Name: WBP2 N-terminal like | |
| Immunogen Catalog Number: AG18347 | Calculated MW: 309 aa, 32 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51240-1, WBP2NL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60839-1-PBS. Detection antibody: 60839-2-PBS. Standard:Ag18347. Range: 0.098-100 ng/mL.