

For Research Use Only

WBP2NL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60839-4-PBS



Basic Information

Catalog Number: 60839-4-PBS	GenBank Accession Number: BC022546	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 164684	CloneNo.: 1F7C12
Source: Mouse	UNIPROT ID: Q6ICG8	
Isotype: IgG1	Full Name: WBP2 N-terminal like	
Immunogen Catalog Number: AG18347	Calculated MW: 309 aa, 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51240-2, WBP2NL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60839-3-PBS. Detection antibody: 60839-4-PBS. Standard:Ag18347. Range: 0.391-100 ng/mL.