

For Research Use Only

TNRC6A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60841-5-PBS



Basic Information

Catalog Number: 60841-5-PBS	GenBank Accession Number: BC024324	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 27327	CloneNo.: 1G3H6
Source: Mouse	UNIPROT ID: Q8NDV7	
Isotype: IgG1	Full Name: trinucleotide repeat containing 6A	
Immunogen Catalog Number: AG34199	Calculated MW: 182 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51246-3, TNRC6A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60841-4-PBS. Detection antibody: 60841-5-PBS. Standard:Ag34199. Range: 1.563-100 ng/mL.