

For Research Use Only

RECK Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60843-5-PBS



Basic Information

Catalog Number:

60843-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26143

GenBank Accession Number:

BC060806

GeneID (NCBI):

8434

UNIPROT ID:

O95980

Full Name:

reversion-inducing-cysteine-rich protein with kazal motifs

Calculated MW:

106 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E7G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51252-3, RECK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60843-4-PBS. Detection antibody: 60843-5-PBS. Standard:Ag26143. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP51252-4, RECK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60843-2-PBS. Detection antibody: 60843-5-PBS. Standard:Ag26143. Range: 0.098-25 ng/mL