

For Research Use Only

Anti-Human STRO-1 (STRO-1)

Catalog Number: 65184-1-Ig



Basic Information

Catalog Number:

65184-1-Ig

Size:

100ug, 0.5 mg/ml

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.



1X10⁶ lysed human blood cells were surface stained with 0.2 ug Anti-Human STRO-1 (65184-1-Ig, Clone:STRO-1) and AlexaFluor-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. CD45 dim, CD34 positive cells were gated.