

For Research Use Only

Anti-Human CD86 (BU63) Rabbit Recombinant Antibody



Catalog Number: 65514-1-RR

Basic Information

Catalog Number: 65514-1-RR	GenBank Accession Number: BC040261	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/mL	GeneID (NCBI): 942	CloneNo.: BU63
Source: Rabbit	ENSEMBL Gene ID: ENSG00000114013	
Isotype: IgG	Full Name: CD86 molecule	
	Calculated MW: 329 aa, 38 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Human

Background Information

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

Storage

Storage:
Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were treated with FC Receptor Block then surface stained with 0.1 ug Anti-Human CD86 Rabbit Recombinant Antibody (65514-1-RR, Clone:BU63) or Rabbit IgG control PolyAb (30000-0-AP), and CoraLite®647-conjugated AffiniPure F(ab')₂ Fragment Donkey Anti-Rabbit IgG (H+L) (SA00014-7). Cells were not fixed. Monocytes were gated.



1x10⁶ human PBMCs were treated with FC Receptor Block then surface stained with 0.1 ug Anti-Human CD86 Rabbit Recombinant Antibody (65514-1-RR, Clone:BU63) or Rabbit IgG control PolyAb (30000-0-AP), and CoraLite®647-conjugated AffiniPure F(ab')₂ Fragment Donkey Anti-Rabbit IgG (H+L) (SA00014-7). Cells were not fixed. Monocytes were gated.