

For Research Use Only

Anti-Rat CD11b (OX42) Mouse IgG2a Recombinant Antibody, PBS Only

Catalog Number: 65606-1-PBS



Basic Information

Catalog Number:

65606-1-PBS

Size:

1mg, 2 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

NM_012711

GeneID (NCBI):

25021

Full Name:

integrin alpha M

Calculated MW:

127 kDa

Purification Method:

Protein A purification

CloneNo.:

OX42

Applications

Tested Applications:

FC

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rat bone marrow cells were surface stained with 0.25 ug Anti-Rat CD11b Mouse IgG2a Recombinant Antibody (65606-1-MR, Clone: OX42) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed. This data was developed using the same antibody clone with 65606-1-PBS in a different storage buffer formulation.



1x10⁶ rat bone marrow cells were surface stained with 0.25 ug Anti-Rat CD11b Mouse IgG2a Recombinant Antibody (65606-1-MR, Clone: OX42) (red) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) (blue) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed. This data was developed using the same antibody clone with 65606-1-PBS in a different storage buffer formulation.