

For Research Use Only

# Anti-Rat CD43 (W3/13) Mouse IgG2a Recombinant Antibody



Catalog Number: 65607-1-MR

## Basic Information

Catalog Number:  
65607-1-MR  
Size:  
100ug, 1 mg/ml  
Source:  
Mouse  
Isotype:  
IgG2a

GenBank Accession Number:  
GeneID (NCBI):  
24796  
Full Name:  
sialophorin  
Calculated MW:  
38 kDa

Purification Method:  
Protein A purification  
CloneNo.:  
W3/13

## Applications

Tested Applications:  
FC  
Species Specificity:  
Rat

## Background Information

### Storage

Storage:  
Store at 2-8°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



$1 \times 10^6$  rat splenocytes were surface stained with 0.25  $\mu$ g Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed.



$1 \times 10^6$  rat splenocytes were surface stained with 0.25  $\mu$ g Anti-Rat CD43 Mouse IgG2a Recombinant Antibody (65607-1-MR, Clone: W3/13) (red) or Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) (blue) and FITC anti-Mouse IgG2a Antibody. Cells were not fixed.