

For Research Use Only

GST Tag Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66001-3-PBS

Basic Information

Catalog Number:

66001-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

FN315687

GeneID (NCBI):

Full Name:

Calculated MW:

26 kDa

Purification Method:

Protein G purification

CloneNo.:

2C6C10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50001-1, GST Tag Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66001-3-PBS. Detection antibody: HRP-conjugated 66001-4-PBS. Standard: gst. Range: 0.156-20ng/mL