

For Research Use Only

GST Tag Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66001-4-PBS

Basic Information

Catalog Number: 66001-4-PBS	GenBank Accession Number: FN315687	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI):	CloneNo.: 1A8C6
Source: Mouse	Full Name:	
Isotype: IgG1	Calculated MW: 26 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50001-1, GST Tag Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66001-3-PBS. Detection antibody: HRP-conjugated 66001-4-PBS. Standard: gst. Range: 0.156-20ng/mL