

For Research Use Only

# HLA class I (HLA-A) Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 66013-3-PBS



## Basic Information

Catalog Number:

66013-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG7488

GenBank Accession Number:

BC003069

GeneID (NCBI):

3105

UNIPROT ID:

P04439

Full Name:

major histocompatibility complex,  
class I, A

Calculated MW:

41 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F8C3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50984-1, HLA class I (HLA-A) Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66013-2-PBS. Detection antibody: 66013-3-PBS. Standard:Ag7488. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP50984-2, HLA class I (HLA-A) Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66013-3-PBS. Detection antibody: 66013-4-PBS. Standard:Ag7488. Range: 3.125-200 ng/mL