

For Research Use Only

MCL1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66026-2-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 66026-2-PBS | GenBank Accession Number: BC107735 | Purification Method: Protein A Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 4170 | CloneNo.: 7D5A3 |
| Source: Mouse | UNIPROT ID: Q07820 | |
| Isotype: IgG2b | Full Name: myeloid cell leukemia sequence 1 (BCL2-related) | |
| Immunogen Catalog Number: AG10609 | Calculated MW: 350 aa, 37 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50655-1, MCL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66026-1-PBS. Detection antibody: 66026-2-PBS. Standard: Ag10609. Range: 0.098-100 ng/mL.