

For Research Use Only

CPT1A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66039-2-PBS

Basic Information

Catalog Number: 66039-2-PBS	GenBank Accession Number: BC000185	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1374	CloneNo.: 6B4A4
Source: Mouse	UNIPROT ID: P50416	
Isotype: IgG3	Full Name: carnitine palmitoyltransferase 1A (liver)	
Immunogen Catalog Number: AG7745	Calculated MW: 88 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50197-1, CPT1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66039-1-PBS. Detection antibody: 66039-2-PBS. Standard: ag7745. Range: 0.195-50 ng/mL