

For Research Use Only

LPCAT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66044-4-PBS



Basic Information

Catalog Number: 66044-4-PBS	GenBank Accession Number: BC020166	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 79888	CloneNo.: 3A5A7
Source: Mouse	UNIPROT ID: Q8NF37	
Isotype: IgG1	Full Name: lysophosphatidylcholine acyltransferase 1	
Immunogen Catalog Number: AG9060	Calculated MW: 534 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50776-4, LPCAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66044-4-PBS. Detection antibody: 66044-2-PBS. Standard: Ag9060. Range: 0.391-100 ng/mL.