

For Research Use Only

NDUFA4L2 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66050-4-PBS



Basic Information

Catalog Number:

66050-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9233

GenBank Accession Number:

BC011910

GeneID (NCBI):

56901

UNIPROT ID:

Q9NRX3

Full Name:

NADH dehydrogenase (ubiquinone) 1
alpha subcomplex, 4-like 2

Calculated MW:

87 aa, 10 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A8F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50529-3, NDUF4L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66050-4-PBS. Detection antibody: 66050-5-PBS. Standard: Ag9233. Range: 0.195-100 ng/mL