

For Research Use Only

NDUFS5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66053-2-PBS



Basic Information

Catalog Number:

66053-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7663

GenBank Accession Number:

BC001884

GeneID (NCBI):

4725

UNIPROT ID:

O43920

Full Name:

NADH dehydrogenase (ubiquinone)

Fe-S protein 5, 15kDa (NADH-coenzyme Q reductase)

Calculated MW:

13 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2D1F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50607-1, NDUFS5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66053-2-PBS. Detection antibody: 66053-1-PBS. Standard: Ag7663. Range: 3.125-100 ng/mL.