

For Research Use Only

CD31 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66065-4-PBS

Basic Information

Catalog Number: 66065-4-PBS	GenBank Accession Number: BC022512	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5175	CloneNo.: 2B1G10
Source: Mouse	ENSEMBL Gene ID: ENSG00000261371	
Isotype: IgG2b	UNIPROT ID: P16284	
Immunogen Catalog Number: AG1787	Full Name: platelet/endothelial cell adhesion molecule	
	Calculated MW: 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50319-1, CD31 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66065-3-PBS. Detection antibody: 66065-4-PBS. Standard: Ag1787. Range: 0.781-200 ng/mL.