

For Research Use Only

Integrin alpha-3 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66070-2-PBS



Basic Information

Catalog Number:

66070-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17293

GenBank Accession Number:

BC150190

GeneID (NCBI):

3675

UNIPROT ID:

P26006

Full Name:

integrin, alpha 3 (antigen CD49C,
alpha 3 subunit of VLA-3 receptor)

Calculated MW:

1066 aa, 119 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3A10G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50516-1, Integrin alpha-3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66070-2-PBS. Detection antibody: 66070-3-PBS. Standard: Ag17293. Range: 0.098-6.25 ng/mL