

For Research Use Only

TELO2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66077-2-PBS



Basic Information

Catalog Number: 66077-2-PBS	GenBank Accession Number: BC017188	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 9894	CloneNo.: 1F12B6
Source: Mouse	UNIPROT ID: Q9Y4R8	
Isotype: IgG1	Full Name: TEL2, telomere maintenance 2, homolog (S. cerevisiae)	
Immunogen Catalog Number: AG8834	Calculated MW: 837 aa, 92 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50016-1, TELO2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66077-2-PBS. Detection antibody: HRP-conjugated 66077-1-PBS. Standard: Ag8834. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP50016-1, TELO2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66077-2-PBS. Detection antibody: 66077-1-PBS. Standard: Ag8834. Range: 1.563-12.5 ng/mL.