

For Research Use Only

# SERPINA3,AACT Monoclonal antibody, PBS Only (Detector)

Catalog Number:66078-3-PBS



## Basic Information

Catalog Number:

66078-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2830

GenBank Accession Number:

BC010530

GeneID (NCBI):

12

UNIPROT ID:

P01011

Full Name:

serpin peptidase inhibitor, clade A  
(alpha-1 antiproteinase, antitrypsin),  
member 3

Calculated MW:

448 aa, 48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E2H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50963-1, SERPINA3,AACT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66078-2-PBS. Detection antibody: 66078-3-PBS. Standard:Ag2830. Range: 3.125-200 ng/mL