

For Research Use Only

# CD22 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66103-2-PBS

## Basic Information

|                                      |                                       |  |
|--------------------------------------|---------------------------------------|--|
| Catalog Number:<br>66103-2-PBS       | GenBank Accession Number:<br>BC109306 | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>933                 | CloneNo.:<br>6D10F4                            |
| Source:<br>Mouse                     | UNIPROT ID:<br>P20273                 |  |
| Isotype:<br>IgG2a                    | Full Name:<br>CD22 molecule           |  |
| Immunogen Catalog Number:<br>AG17986 | Calculated MW:<br>847 aa, 95 kDa      |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50311-1, CD22 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66103-2-PBS. Detection antibody: 66103-3-PBS. Standard: Ag17986. Range: 0.781-100 ng/mL.