

For Research Use Only

IL-16 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66145-2-PBS



Basic Information

Catalog Number:

66145-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG16466

GenBank Accession Number:

BC040272

GeneID (NCBI):

3603

UNIPROT ID:

Q14005

Full Name:

interleukin 16 (lymphocyte chemoattractant factor)

Calculated MW:

130 aa, 15 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3D12C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51113-1, IL-16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66145-1-PBS. Detection antibody: 66145-2-PBS. Standard: Ag16466. Range: 3.125-200 ng/mL.