

For Research Use Only

IL6 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66146-2-PBS

Basic Information

Catalog Number: 66146-2-PBS	GenBank Accession Number: BC015511	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3569	CloneNo.: 3B11C11
Source: Mouse	ENSEMBL Gene ID: ENSG00000136244	
Isotype: IgG2b	UNIPROT ID: P05231	
Immunogen Catalog Number: AG10646	Full Name: interleukin 6 (interferon, beta 2)	
	Calculated MW: 212 aa, 24 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50032-2, IL-6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68659-3-PBS. Detection antibody: 66146-2-PBS. Standard:HZ-1019. Range: 1.563-100 ng/mL.