

For Research Use Only

Caspase 9 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66169-2-PBS

Basic Information

Catalog Number: 66169-2-PBS	GenBank Accession Number: BC002452	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 842	CloneNo.: 1C1H12
Source: Mouse	UNIPROT ID: P55211	
Isotype: IgG1	Full Name: caspase 9, apoptosis-related cysteine peptidase	
Immunogen Catalog Number: AG20813	Calculated MW: 46 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50283-1, Caspase 9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66169-1-PBS. Detection antibody: 66169-2-PBS. Standard: Ag20813. Range: 1.563-200 ng/mL.