

For Research Use Only

# CD71 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66180-2-PBS



## Basic Information

<b>Catalog Number:</b> 66180-2-PBS	<b>GenBank Accession Number:</b> BC001188	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 7037	<b>CloneNo.:</b> 5F6H4
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000072274	
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> P02786	
<b>Immunogen Catalog Number:</b> AG21612	<b>Full Name:</b> transferrin receptor (p90, CD71)	
	<b>Calculated MW:</b> 85 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50414-1, CD71 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68962-1-PBS. Detection antibody: 66180-2-PBS. Standard:Eg0339. Range: 1.563-50 ng/mL