

For Research Use Only

KRT18 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66187-2-PBS

Basic Information

Catalog Number:

66187-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1260

GenBank Accession Number:

BC000180

GeneID (NCBI):

3875

UNIPROT ID:

P05783

Full Name:

keratin 18

Calculated MW:

48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

5D10E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50238-1, KRT18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66187-2-PBS. Detection antibody: 66187-1-PBS. Standard: Ag1260. Range: 0.391-100 ng/mL.