

For Research Use Only

ERK1/2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66192-2-PBS



Basic Information

Catalog Number: 66192-2-PBS	GenBank Accession Number: BC013992	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5595	CloneNo.: 3C2E12
Source: Mouse	UNIPROT ID: P27361	
Isotype: IgG1	Full Name: mitogen-activated protein kinase 3	
Immunogen Catalog Number: AG20350	Calculated MW: 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50383-1, ERK1/2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67170-2-PBS. Detection antibody: 66192-2-PBS. Standard:Ag20350. Range: 3.125-50 ng/mL



Cytometric bead array standard curve of MP50383-2, ERK1/2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66192-3-PBS. Detection antibody: 66192-2-PBS. Standard:Ag20350. Range: 1.563-25 ng/mL