

For Research Use Only

IDH1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66197-2-PBS



Basic Information

Catalog Number: 66197-2-PBS	GenBank Accession Number: BC012846	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3417	CloneNo.: 4B11F2
Source: Mouse	UNIPROT ID: O75874	
Isotype: IgG1	Full Name: isocitrate dehydrogenase 1 (NADP+), soluble	
Immunogen Catalog Number: AG19293	Calculated MW: 414 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51184-1, IDH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66197-2-PBS. Detection antibody: 66197-1-PBS. Standard: Ag19293. Range: 3.125-200 ng/mL.