

For Research Use Only

# NINJ2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66202-3-PBS



## Basic Information

Catalog Number: 66202-3-PBS	GenBank Accession Number: BC057766	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4815	CloneNo.: 1D1B10
Source: Mouse	UNIPROT ID: Q9NZG7	
Isotype: IgG2b	Full Name: ninjurin 2	
Immunogen Catalog Number: AG5343	Calculated MW: 142 aa, 16 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50693-1, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-2-PBS. Detection antibody: 66202-3-PBS. Standard:Ag5343. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50693-2, NINJ2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66202-4-PBS. Detection antibody: 66202-3-PBS. Standard:Ag5343. Range: 0.391-100 ng/mL.