

For Research Use Only

# MYOD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66214-3-PBS



## Basic Information

Catalog Number: 66214-3-PBS	GenBank Accession Number: BC064493	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4654	CloneNo.: 3F6F5
Source: Mouse	UNIPROT ID: P15172	
Isotype: IgG1	Full Name: myogenic differentiation 1	
Immunogen Catalog Number: AG13512	Calculated MW: 320 aa, 35 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51167-1, MYOD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66214-2-PBS. Detection antibody: 66214-3-PBS. Standard: Ag13512. Range: 12.5-200 ng/mL