

For Research Use Only

CTNNA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66221-2-PBS



Basic Information

Catalog Number: 66221-2-PBS	GenBank Accession Number: BC031262	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1495	CloneNo.: 1F8B7
Source: Mouse	UNIPROT ID: P35221	
Isotype: IgG2b	Full Name: catenin (cadherin-associated protein), alpha 1, 102kDa	
Immunogen Catalog Number: AG23603	Calculated MW: 906 aa, 100 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51132-1, CTNNA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66221-2-PBS. Detection antibody: 66221-3-PBS. Standard: Ag23603. Range: 0.781-100 ng/mL.