

For Research Use Only

SALL4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 66236-2-PBS

Basic Information

Catalog Number:

66236-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17480

GenBank Accession Number:

BC111714

GeneID (NCBI):

57167

UNIPROT ID:

Q9UJQ4

Full Name:

sal-like 4 (Drosophila)

Calculated MW:

1053 aa, 112 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50336-1, SALL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66236-2-PBS. Detection antibody: 66236-3-PBS. Standard: Ag17480. Range: 1.563-100 ng/mL.