

For Research Use Only

SALL4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66236-3-PBS

Basic Information

Catalog Number: 66236-3-PBS	GenBank Accession Number: BC111714	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 57167	CloneNo.: 4D5E12
Source: Mouse	UNIPROT ID: Q9UJQ4	
Isotype: IgG1	Full Name: saI-like 4 (Drosophila)	
Immunogen Catalog Number: AG17480	Calculated MW: 1053 aa, 112 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50336-1, SALL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66236-2-PBS. Detection antibody: 66236-3-PBS. Standard:Ag17480. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50336-2, SALL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66236-4-PBS. Detection antibody: 66236-3-PBS. Standard:Ag17480. Range: 1.563-100 ng/mL.