

For Research Use Only

CA9 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66243-3-PBS

Basic Information

Catalog Number:

66243-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1540

GenBank Accession Number:

BC014950

GeneID (NCBI):

768

UNIPROT ID:

Q16790

Full Name:

carbonic anhydrase IX

Calculated MW:

459 aa, 50 kDa

Purification Method:

Protein G purification

CloneNo.:

1B12G2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50143-1, CA9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66243-2-PBS. Detection antibody: 66243-3-PBS. Standard: Ag1540. Range: 1.563-100 ng/mL.