

For Research Use Only

CA9 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66243-5-PBS

Basic Information

Catalog Number: 66243-5-PBS	GenBank Accession Number: BC014950	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 768	CloneNo.: 1H9F10
Source: Mouse	UNIPROT ID: Q16790	
Isotype: IgG1	Full Name: carbonic anhydrase IX	
Immunogen Catalog Number: AG1540	Calculated MW: 459 aa, 50 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50143-2, CA9 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 66243-4-PBS. Detection antibody: 66243-5-PBS. Standard: Ag1540. Range: 1.563-100 ng/mL