

For Research Use Only

C-Reactive Protein/CRP Monoclonal antibody

Catalog Number: 66250-1-Ig

Featured Product

12 Publications



Basic Information

Catalog Number:

66250-1-Ig

Concentration:

3200 ug/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9883

GenBank Accession Number:

BC020766

GeneID (NCBI):

1401

UNIPROT ID:

P02741

Full Name:

C-reactive protein, pentraxin-related

Calculated MW:

224 aa, 25 kDa

Observed MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

3F10A11

Recommended Dilutions:

WB: 1:2000-1:10000

IHC: 1:50-1:500

IF-P: 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF-P, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat, pig

Cited Species:

human, mouse, rat, monkey

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : serum from mouse injected with bacteria, pig liver tissue, rat liver tissue, serum from mouse injected with bacteria tissue

IHC : human liver tissue,

IF-P : human liver tissue,

Background Information

C-reactive protein (CRP) is an acute-phase serum protein synthesized by the liver in response to interleukin-6 (IL-6) during inflammation. The name of CRP derives from its ability to react with the C polysaccharide of *Streptococcus pneumoniae*. CRP is an annular, pentameric protein that belongs to the pentraxin family of proteins. CRP displays several functions associated with host defense: it promotes agglutination, bacterial capsular swelling, phagocytosis and complement fixation through its calcium-dependent binding to phosphorylcholine. It is used mainly as a marker of inflammation.

Notable Publications

Author	Pubmed ID	Journal	Application
Amrita K Cheema	35675853	Int J Radiat Oncol Biol Phys	WB
Shaoning Li	33599398	Brain Behav	WB, IHC
Anika Alberts	33363539	Front Immunol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Serum from mouse injected with bacteria were subjected to SDS PAGE followed by western blot with 66250-1-Ig (C-Reactive Protein/CRP Antibody) at dilution of 1:2000. Horseradish peroxidase conjugated protein A was used for signal developing instead of secondary antibody incubated at room temperature for 1.5 hours.




Serum from normal mouse (left) and mouse injected with bacteria (right) were subjected to SDS PAGE followed by western blot with 66250-1-Ig (C-Reactive Protein/CRP Antibody) at dilution of 1:1000. Horseradish peroxidase conjugated protein A was used for signal developing instead of secondary antibody incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 66250-1-Ig (C-Reactive Protein/CRP Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 66250-1-Ig (C-Reactive Protein/CRP Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (4% PFA) fixed human liver tissue using 66250-1-Ig (CRP antibody) at dilution of 1:100 and CoraLite488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).