

For Research Use Only

PHOX2B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66254-3-PBS



Basic Information

Catalog Number:

66254-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13196

GenBank Accession Number:

BC017199

GeneID (NCBI):

8929

UNIPROT ID:

Q99453

Full Name:

paired-like homeobox 2b

Calculated MW:

31 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E12C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP51004-1, PHOX2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66254-2-PBS. Detection antibody: 66254-3-PBS. Standard:Ag13196. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51004-3, PHOX2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66254-6-PBS. Detection antibody: 66254-3-PBS. Standard:Ag13196. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51004-4, PHOX2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66254-5-PBS. Detection antibody: 66254-3-PBS. Standard:Ag13196. Range: 0.391-100 ng/mL