

For Research Use Only

VCAM-1/CD106 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66294-2-PBS



Basic Information

Catalog Number: 66294-2-PBS	GenBank Accession Number: BC017276	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7412	CloneNo.: 1G12C5
Source: Mouse	UNIPROT ID: P19320	
Isotype: IgG1	Full Name: vascular cell adhesion molecule 1	
Immunogen Catalog Number: AG24697	Calculated MW: 739 aa, 81 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50328-1, VCAM-1/CD106 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66294-2-PBS. Detection antibody: 66294-3-PBS. Standard: Ag24697. Range: 0.781-12.5 ng/mL