

For Research Use Only

EPCAM/CD326 Monoclonal antibody, PBS Only (Capture)

Catalog Number:66316-6-PBS



Basic Information

Catalog Number: 66316-6-PBS	GenBank Accession Number: BC014785	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4072	CloneNo.: 3G10C12
Source: Mouse	UNIPROT ID: P16422	
Isotype: IgG2a	Full Name: epithelial cell adhesion molecule	
Immunogen Catalog Number: AG15393	Calculated MW: 314 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50309-4, EPCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66316-6-PBS. Detection antibody: 66316-7-PBS. Standard: Ag15393. Range: 3.125-100 ng/mL.