

For Research Use Only

Peripherin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66317-3-PBS



Basic Information

Catalog Number:

66317-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG11116

GenBank Accession Number:

BC032703

GeneID (NCBI):

5630

UNIPROT ID:

P41219

Full Name:

peripherin

Calculated MW:

470 aa, 54 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B7C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51128-1, Peripherin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66317-2-PBS. Detection antibody: 66317-3-PBS. Standard:Ag11116. Range: 1.563-100 ng/mL