

For Research Use Only

SRD5A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 66329-2-PBS



Basic Information

Catalog Number: 66329-2-PBS	GenBank Accession Number: BC007033	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6715	CloneNo.: 1B5H10
Source: Mouse	UNIPROT ID: P18405	
Isotype: IgG2b	Full Name: steroid-5-alpha-reductase, alpha polypeptide 1 (3-oxo-5 alpha-steroid delta 4-dehydrogenase alpha 1)	
Immunogen Catalog Number: AG23351		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50930-1, SRD5A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66329-1-PBS. Detection antibody: 66329-2-PBS. Standard: Ag23351. Range: 0.781-100 ng/mL.