

For Research Use Only

# BCL6 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 66340-2-PBS

## Basic Information

Catalog Number:

66340-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG15519

GenBank Accession Number:

BC150184

GeneID (NCBI):

604

UNIPROT ID:

P41182

Full Name:

B-cell CLL/lymphoma 6

Calculated MW:

706 aa, 79 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H5B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50215-1, BCL6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66340-5-PBS. Detection antibody: 66340-2-PBS. Standard: Ag15519. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP80008-1, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-1-PBS. Detection antibody: 81751-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL.