

For Research Use Only

# MKP-2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 66349-3-PBS

## Basic Information

Catalog Number:

66349-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1206

GenBank Accession Number:

BC014565

GeneID (NCBI):

1846

UNIPROT ID:

Q13115

Full Name:

dual specificity phosphatase 4

Calculated MW:

43 kDa

Purification Method:

Protein G purification

CloneNo.:

1D10C12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50101-1, MKP-2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66349-2-PBS. Detection antibody: 66349-3-PBS. Standard: Ag1206. Range: 6.25-100 ng/mL.