

For Research Use Only

MMP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66366-2-PBS



Basic Information

Catalog Number: 66366-2-PBS	GenBank Accession Number: BC002576	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4313	CloneNo.: 2A8A11
Source: Mouse	UNIPROT ID: P08253	
Isotype: IgG1	Full Name: matrix metalloproteinase 2 (gelatinase A, 72kDa gelatinase, 72kDa type IV collagenase)	
Immunogen Catalog Number: AG25039	Calculated MW: 72 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51054-1, MMP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66366-2-PBS. Detection antibody: 66366-3-PBS. Standard:Ag25039. Range: 0.098-100 ng/mL.