

For Research Use Only

# MMP2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66366-4-PBS



## Basic Information

Catalog Number:

66366-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25039

GenBank Accession Number:

BC002576

GeneID (NCBI):

4313

UNIPROT ID:

P08253

Full Name:

matrix metalloproteinase 2  
(gelatinase A, 72kDa gelatinase,  
72kDa type IV collagenase)

Calculated MW:

72 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H3F11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51054-2, MMP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66366-4-PBS. Detection antibody: 66366-1-PBS. Standard:Ag25039. Range: 0.098-100 ng/mL.