

For Research Use Only

# TrkC Monoclonal antibody, PBS Only

Catalog Number: 66380-1-PBS



## Basic Information

Catalog Number:

66380-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2635

GenBank Accession Number:

BC013693

GeneID (NCBI):

4916

UNIPROT ID:

Q16288

Full Name:

neurotrophic tyrosine kinase,  
receptor, type 3

Calculated MW:

839 aa, 94 kDa

Observed MW:

110 kDa, 160 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7F11

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat, mouse, pig

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66380-1-Ig (TrkC) antibody at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66380-1-PBS in a different storage buffer formulation.