

For Research Use Only

OGN Monoclonal antibody, PBS Only (Capture)

Catalog Number: 66382-3-PBS



Basic Information

Catalog Number:

66382-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3484

GenBank Accession Number:

BC037273

GeneID (NCBI):

4969

UNIPROT ID:

P20774

Full Name:

osteoglycin

Calculated MW:

298 aa, 34 kDa

Purification Method:

Protein G purification

CloneNo.:

1E4E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50723-2, OGN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66382-3-PBS. Detection antibody: 66382-4-PBS. Standard: Ag3484. Range: 0.098-100 ng/mL.